

**Annex 2:** description of the kits used at the CIRCB laboratory

<b>Kit name</b>	<b>Manufacturer</b>	<b>Method</b>	<b>Prequalified</b>	<b>Sensitivity IgG/IgM % (CI%)</b>	<b>Specificity IgG/IgM % (CI %)</b>
TOXO IgG/IgM Rapid Test Cassette	JusChek	ICT	IND	96.0/94,0 (86.3-99.5)/ (83.5-98.7)	98.5/98.8 (98.6-99.4) / (97.1-99.6)
Toxoplasma IgG/IgM ELISA	DIALAB	EIA	IND	> 99.9 / > 99.9 (93,4-100) / (88.1-100)	99.0/99.0 (98.2-99.4)/ (98.3-99.5)
VIDAS®TOXO IgG II (TXG) / IgM (TXM)	Biomérieux	ELFA	IND	99,65/96,00 (98.75-99.96)/ (91.43-98.18)	99.92/99.25 (99.58-100)/ (98.66-99.59)

CI= Confidence interval, CIRCB=Chantal Biya International Reference Center, EIA= Enzyme Immunoassay, ELFA= Enzyme-linked fluorescent assay, ELISA=Enzyme linked immune sorbent assay, ICT= Immunochroma-tographic, Ig=Immunoglobulin, IND=Indeterminate, TX=Toxo